

Accent *Prime*
Love Your Shape

Accentuate

Hands-free Treatment
The Freedom to Do More

Alma

To achieve optimal results, the recommended treatment protocol for skin tightening and body contouring is a combination of both UltraSpeed and Accentuate applicators. UltraSpeed uses a new type of Ultrasound technology delivered with an extra-large applicator plate for high speed fat reduction. This applicator features a plate type sonotrode which emits simultaneously longitudinal and transverse ultrasonic waves through concentric profiles.

The US energy is distributed homogenously throughout the treatment area, effectively disrupting stubborn fat cells. Following the use of UltraSpeed, Accentuate complements the treatment by delivering RF energy into defined areas, heating tissues at varying depths, and causing a thermal effect within the skin. Accentuate contours and shapes the treated area by tightening the skin after the fat cells have already been disrupted by UltraSpeed. This synergetic pairing of technologies and applicators, addressing both fat reduction and skin tightening in a single treatment session, offers the ideal solution for body contouring.

A treatment protocol composed of 20 minutes of UltraSpeed treatment followed by 40 minutes of Accentuate, provides highly effective results for large areas.



Accentuate







UltraSpeed

APPLICATOR HIGHLIGHTS

- ✓ AlmaWave 40.68 MHz
- ✓ Powerful - up to 300 Watts
- ✓ Real-time patient temperature control
- ✓ Ergonomic design
- ✓ Remote heating
- ✓ Independent single or double applicators use
- ✓ A washable mattress and additional hygienic covers

TREATMENT AREAS

-  Abdomen
-  Outer thighs
-  Inner thighs
-  Buttocks



Accentuate offers a significant advantage for practitioners given the fact it is a hands-free treatment, making it exceptionally time-efficient and valuable. Practitioners can position the applicator on the patient, define the treatment parameters and benefit from the ability to conduct multiple treatments at once. Accentuate distributes RF-energy homogenously across the treated area in close proximity to the skin, yet, without direct contact, thus ensuring efficacy while maintaining treatment safety.



Throughout the treatment an internal cooling fan is continuously operated to ensure that patient's skin remains dry and comfortable at all times. Temperature is continuously monitored and stabilized, making sure that the treatment is safe even when the practitioner is not present. Accentuate is a flexible device and can be deployed using one or two of its applicators for consistent coverage of large treatment areas. In addition, it features built-in pre-defined parameters to suit patients' individual needs.



THE ALMAWAVE EFFECT

The AlmaWave provides deep, homogenous dielectric heating for visible, long lasting results. This intensive RF energy induces heating of various layers of the skin, contracting collagen fibers and stimulating the formation of new collagen to improve skin's thickness and alignment. Utilizing the optimal AlmaWave 40.68MHz wavelength enables Accentuate to deliver energy deeper and more effectively within the skin's tissue.

Thermal Imaging



HANDS-FREE FOR ULTIMATE EFFICIENCY

In today's world, the ability to manage resources and multi-task plays a huge role when it comes to maintaining desired levels of efficiency and meeting ongoing business goals. This is particularly true for technology-driven practitioners, especially when they want to develop their practice and offer top-notch customer service. The chance to conduct multiple treatments simultaneously, while providing excellent care is undoubtedly the key to professional leadership and business growth.

Accentuate is a hands-free, stationary applicator designed for body contouring and skin tightening, incorporating resonant distributed RF energy and exclusively used with Alma's Accent Prime platform. Utilizing the AlmaWave 40.68 MHZ technology, Accentuate shapes and tightens large areas of the body safely and effectively by delivering energy within the skin.

Offering independent left & right side parameters, Accentuate also enables real-time patient temperature monitoring, adding increased control and freedom for practitioners and patients alike. Accentuate provides optimal results when used in a combined protocol following UltraSpeed. This synergetic combination of RF and Ultrasound technologies enables enhanced effectiveness and significantly improved results.



BENEFITS

-  HANDS-FREE TREATMENT
-  COMBINED TECHNOLOGIES FOR ENHANCED RESULTS
-  SAFE AND EFFECTIVE
-  CLINICALLY PROVEN
-  NO DOWNTIME; NO SIDE EFFECTS
-  COMFORTABLE, PAINLESS TREATMENT
-  HOMOGENEOUS COVERAGE OF LARGE AREAS
-  ENHANCED BUSINESS PRODUCTIVITY

“I find the protocol that combines both UltraSpeed and Accentuate leads to the best treatment results for large areas of the body – including the abdomen, thighs and buttocks. In just 1 hour, my patients benefit from tighter and contoured skin. The added value with Accentuate is my ability to treat several patients at once without having to be present personally during the treatment.”

Dr. Stefania Enginoli, Studio Medico, Milano



Accent *Prime*
Love Your Shape

Accentuate

Hands-free Treatment
The Freedom to Do More

Alma Lasers GmbH

Nordostpark 100-102
90411 Nürnberg, Germany
Tel. + 49-911 / 89-11-29-0
Fax + 49-911 / 89-11-29-99
info@alma-lasers.de
www.almalasers.com

© 2018 Alma Lasers.

All rights reserved. Alma Lasers, its logo, Alma ,Accent Prime, UltraSpeed, Accentuate, AlmaWave, Unipolar are trademarks and/or registered trademarks of Alma Lasers in the United States and/or other countries.

info@almalasers.com
www.almaaccentprime.com

Connect to Alma



PBAC22021801_01

